

Extraction of Cocaine and Metabolites From Urine Using ISOLUTE SLE+ prior to LC-MS/MS Analysis

Introduction

This application note describes the extraction of cocaine and a full range of cocaine metabolites from human urine prior to analysis and quantitation by LC-MS/MS




Figure 1. Structure of Benzoylecgonine

ISOLUTE SLE+ Supported Liquid Extraction plates and columns offer an efficient alternative to traditional liquid-liquid extraction (LLE) for bioanalytical sample preparation, providing high analyte recoveries, no emulsion formation, and significantly reduced sample preparation

Analytes

Cocaine, Norcocaine, Cocaethylene, Benzoylecgonine, Ecgonine Methyl Ester (EME) and Anhydroecgonine Methyl Ester (AEME).

Sample Preparation Procedure

Sample pre-treatment: Dilute urine (1 mL) with 0.1% ammonium hydroxide (aq) (1 mL). Spike ISTD (10µL of BZE-d3; total concentration 50ng/mL). Vortex mix thoroughly.

ISOLUTE SLE+ 200 µL Supported Liquid Extraction Plate, part number 820-0200-P01

Sample loading: Load the pre-treated urine (200 µL) onto the plate and apply a pulse of vacuum or positive pressure to initiate flow. Allow the sample to adsorb for 5 minutes.

Analyte extraction: Apply DCM/IPA (95:5, v/v) (1 mL) and allow to flow under gravity for 5 minutes. Apply vacuum or positive pressure to pull through any remaining extraction solvent.

ISOLUTE SLE+ 1 mL Sample Volume columns, part number 820-0140-C

Sample loading: Load the pre-treated urine (1 mL) onto the column and apply a pulse of vacuum or positive pressure to initiate flow. Allow the sample to adsorb for 5 minutes.

Analyte extraction: Apply DCM/IPA (95:5, v/v) (2.5 mL) and allow to flow under gravity for 5 minutes. Apply a further aliquot of DCM/IPA (95:5, v/v) (2.5 mL) and allow to flow for another 5 minutes. Apply vacuum or positive pressure to pull through any remaining extraction solvent.

Post-extraction: Add 50mM HCl in MeOH (100 µL) to stabilize AEME and EME during evaporation. Evaporate to dryness at 40 °C.

Reconstitution: **820-0200-P01:** Water/Methanol (90/10, v/v) (200 µL)

820-0140-C: Water/Methanol (90/10, v/v) (500 µL)

UPLC Conditions

Instrument: Waters Acquity UPLC fitted with 20 µL loop

Column: Waters UPLC BEH C18 (1.7 µm, 100 x 2.1 mm id)

Sample temperature: 20 °C

Column temperature: 40 °C

Injection volume: 10 µL (partial loop with overfill)

Mobile phase: A= 0.1% ammonium hydroxide (aq)
B= 0.1% ammonium hydroxide/MeOH

Gradient:

Time (min)	%A	%B	Curve
0	60	40	-
2	10	90	6
3	60	40	11

Mass Spectrometry Conditions

Instrument: Premier XE triple quadrupole mass spectrometer equipped with an electrospray interface for mass analysis.

Desolvation temperature: 450 °C

Ion source temperature: 150 °C

Table 1. Positive ions acquired in the multiple reaction monitoring (MRM) mode

Analyte	Transition	Dwell Time (s)	Cone voltage (V)	Collision energy (eV)
EME	200 > 81.9 QUANT	0.07	30	23
	200 > 182.0 QUAL			17
BZE	290.1 > 168.0 QUANT	0.07	30	18
	290.1 > 149.9 QUAL			22
BZE-d3 ISTD	293.1 > 171.0	0.07	30	18
AEME	182.0 > 122.0 QUANT	0.2	30	19
	182.0 > 90.9 QUAL			21
Norcocaine	290.1 > 135.9 QUANT	0.05	30	22
	290.1 > 107.9 QUAL			26
Cocaine	304.2 > 182.0 QUANT	0.05	30	20
	304.2 > 149.9 QUAL			25
Cocaethylene	318.1 > 196.0 QUANT	0.05	30	19
	318.1 > 149.9 QUAL			23

Results

This ISOLUTE SLE+ protocol demonstrates analyte recoveries ranges from 67-100% as shown in figure 2. RSDs were below 5% for all analytes. Figure 3 shows an example calibration curve for the cocaine metabolite, benzoylecgonine. Figure 4 shows the MRM transitions, both quantification and qualification ions, for analytes extracted using the supported liquid extraction protocol.




Figure 2. Typical analyte % recoveries for a extracted Cocaine and metabolites (n=7) using the ISOLUTE SLE+ protocol

Compound name: BZE Quant
 Correlation coefficient: $r = 0.999740$, $r^2 = 0.999479$
 Calibration curve: $1.07691 * x + 0.0364176$
 Response type: Internal Std (Ref 5), Area * (IS Conc. / IS Area)
 Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None

Compound name: BZE Qual
 Correlation coefficient: $r = 0.999554$, $r^2 = 0.999108$
 Calibration curve: $0.0988376 * x + 0.0152144$
 Response type: Internal Std (Ref 5), Area * (IS Conc. / IS Area)
 Curve type: Linear, Origin: Exclude, Weighting: 1/x, Axis trans: None




Figure 3. Example calibration curve for the cocaine metabolite Benzoylcegonine

Table 2. Estimated Limits of Quantitation for all analytes using the ISOLUTE SLE+ 1 mL sample volume column

Analyte	Estimated LOQ
EME	< 100 pg/mL
BZE	< 100 pg/mL
AEME	< 200 pg/mL
Norcocaine	< 400 pg/mL, restricted by qualifier ion
Cocaine	< 25 pg/mL
Cocaethylene	< 25 pg/mL




Figure 4. MRM transitions chromatograms for Cocaine and metabolites at 1 ng/mL

Ordering information

Part number	Description	Quantity
820-0200-P01	ISOLUTE SLE+ 200 Supported Liquid Extraction Plate	1
820-0140-C	ISOLUTE SLE+ 1 mL Sample Volume columns	30
820-0140-CG	ISOLUTE SLE+ 1 mL Sample Volume columns (tablets for use with PPM-48)	30
PPM-96	PRESSURE+96 Positive Pressure Manifold 96 Well	1
PPM-48	PRESSURE+48 Positive Pressure Manifold 48 position	1
121-9600	VacMaster 96 Sample Processing Manifold	1
121-1016	VacMaster-10 Sample Processing Manifold	1
121-2016	VacMaster-20 Sample Processing Manifold	1
SD-9600-DHS-UK	SPE Dry 96 240 V	1

To search and download more of Biotage's extensive database of sample preparation applications please visit <http://www.biotage.com/applications> or scan the QR code with your smart phone to go direct.



NORTH AMERICA

Main Office: +1 704 654 4900
Toll Free: +1 800 446 4752
Fax: +1 704 654 4917
Order Tel: +1 704 654 4900
Order Fax: +1 434 296 8217
ordermailbox@biotage.com
US-1-pointsupport@biotage.com

EUROPE

Main Office: +46 18 56 5900
Fax: +46 18 59 1922
Order Tel: +46 18 56 57 10
Order Fax: +46 18 56 57 05
order@biotage.com
EU-1-pointsupport@biotage.com

JAPAN

Tel: +81 3 5627 3123
Fax: +81 3 5627 3121
jp_order@biotage.com
JP-1-pointsupport@biotage.com

CHINA

Tel: +86 21 2898 6655
Fax: +86 21 2898 6153
CN_order@biotage.com
CN-1-pointsupport@biotage.com
www.biotage.com

Distributors

To locate a distributor please visit our Web site at www.biotage.com